

Room SafeGuard

Cat.No.: P-CC-001

Size: 480mL

Product Introduction

Product Overview In the process of cell culture, contamination is a headache for researchers. Among the contamination sources, researchers often ignore bacteria and fungi in cell rooms, incubators, water baths, and water trays. Based on accumulated experience for years in the field of cell culture, Pricella has independently developed a set of solutions specific to environmental contamination, including cell room contamination, incubator contamination, water bath contamination, water tray contamination, etc. With a strong function of inactivating bacteria and mold spores, Room SafeGuard is mainly used for the disinfection of fungi and bacteria in the cell culture room and CO2 incubators. It is colorless and has a clear fragrance, and it can be sprayed directly without any toxic effects on the human body. By removing bacteria, fungi, and their spores in the cell room space and the incubator, it can effectively prevent microbial contamination, thereby reducing the waste of effort, money, time, and experimental materials caused by cell contamination.

Matters Need Attention

1. After disinfection on cell culture room and cell incubator with Room SafeGuard, seal the disinfected space for 15-30 minutes before reuse. It is suggested to use it with the daily ultraviolet radiation for better effect.; 2. If the Room SafeGuard contacts with the skin, eyes, etc., wash immediately with plenty of water, and seek medical treatment immediately if necessary.; 3. Wear a lab coat and disposable gloves for operation.

Specifications

Volume	480mL
Form	Liquid
Animal Origin Ingredient	Without
Valid Period	24 months
Storage Condition	2~8°C, shading light
Shipping Condition	RT